



Knee Extension Assist. Model 9000 GX-Knee

Based on the design of a traditional high-quality Becker knee joint, the **GX-Knee** comes equipped with a pneumatic spring on the lateral joint to assist with knee extension during the swing phase of gait. A medial knee joint is also included to help control knee alignment.

The **GX-Knee** is a free motion knee joint that is designed to provide extension assist, not prevent knee flexion.



9000

Features:

- Available with your choice of a 125N pneumatic spring (for standard applications), or 175N (for high activity patients).
- Free knee with pneumatic extension assist.
- Paired with a medial joint to control knee alignment.
- Available with aluminum, stainless steel or titanium uprights.

**Becker Orthopedic –
always a step ahead.**



Lateral joint with
pneumatic spring



(800) 521-2192 (248) 588-7480

BeckerOrthopedic.com