



## Security. Freedom. Convenience. Model 2009 Automatic Ring Lock Knee Joint

**Model 2009** is a non-protrusion knee joint that combines the simplicity and security of a standard drop lock with automatic locking.

When the joints reach full extension the rings are automatically engaged by springs that are internally housed within the midsections. This feature gives the user the convenience of automatic locking.

For 75 years Becker Orthopedic has led the industry with high quality,




2009

innovative orthotic components and **central fabrication services.**

This tradition continues today. Whether you are looking for conventional orthotic knee joints (**Models 1002** and **1003**), or the latest **Model 2009**, Becker Orthopedic has the right orthotic option for you.

**Becker Orthopedic –  
always a step ahead.**



1002

1003



©2008 Becker Orthopedic Appliance Co. All rights reserved.



(800) 521-2192  
(248) 588-7480  
[BeckerOrthopedic.com](http://BeckerOrthopedic.com)