



## Easily Adjust A-P Alignment. Model 1013 Posterior Offset Adjustable Knee Joint

**Model 1013** is a posterior offset adjustable knee joint that allows the orthotist to fine tune the A-P alignment of the knee; these features enhance knee stability and help broaden the clinical application.

The A-P alignment and the extension/flexion angle can be changed by simply removing two screws; this allows easy adjustability without the need to physically alter the angle of the uprights.



1013

The **1013** comes equipped with an optional bail lock that can be engaged when the patient reaches full extension.

Depending on patient characteristics and biomechanical requirements, the locking levers can be removed and replaced with the included plugs, allowing the joint to function as a free motion posterior offset knee joint.

Becker Orthopedic -  
always a step ahead.



1013 with  
locking lever

 **BECKER**  
*A Tradition of Excellence A Commitment to Innovation*

(800) 521-2192 (248) 588-7480

[BeckerOrthopedic.com](http://BeckerOrthopedic.com)