

# safetySTRIDE™

Natural Motion. Mechanical Stability.



**Natural Motion:** The SafetyStride™ automatically locks and unlocks during the gait cycle, providing natural movement with superior stability and comprehensive knee joint control.

**Mechanical Stability:** Designed to unlock at terminal stance, an internal lever re-engages during swing phase to ensure knee joint stability prior to heel contact. Individuals who intermittently fail to reach full extension will now have the added security and stability they require while ambulating.

The SafetyStride™ is fully interchangeable with other joints in the Becker Stride Control family.



The SafetyStride™ orthotic knee joint should be considered for individuals with significant quadriceps deficit and generalized lower extremity weakness.

US Patent 7,462,159 B1

## safetySTRIDE™

natural motion. mechanical stability.

### FullStride™ & SafetyStride™ with NEW GX-Assist Option

The FullStride™ and SafetyStride™ may now also be ordered with a GX-Assist option that utilizes a pneumatic spring to assist in knee extension during the swing phase of gait.

For details and ordering information, please contact our customer service department and ask for a copy of our new **Stance Control Overview Guide**.

 **BECKER**  
*A Tradition of Excellence A Commitment to Innovation*

(800) 521-2192 (248) 588-7480  
[BeckerOrthopedic.com](http://BeckerOrthopedic.com)