



## Secure. Convenient. Modular Design. Model 1017 Automatic Angled Lever Lock Knee Joint

**Model 1017** is a modular lever lock knee joint with a locking mechanism that will automatically engage when the joint reaches full extension.

The long angled lever is easily operated by the patient and can be shortened or contoured to suit individual needs. In addition, **Model 1017** comes with an optional anterior pull-cable that allows the patient to simultaneously disengage both locks when necessary.




1017

For over 75 years, Becker Orthopedic has led the industry with high quality, innovative orthotic components and **central fabrication services**.

This tradition continues today. Whether you are looking for conventional orthotic knee joints (**Models 1002 and 1003**), or the latest **Model 1017**, Becker Orthopedic has the right orthotic option for you.

**Becker Orthopedic –  
always a step ahead.**



1002

1003

**BECKER**  
A Tradition of Excellence A Commitment to Innovation

(800) 521-2192 (248) 588-7480

[BeckerOrthopedic.com](http://BeckerOrthopedic.com)

ISO 9001,  
CERTIFIED



©2009 Becker Orthopedic Appliance Co. All rights reserved.